

Hearing Aid Care Leaflet

Date of Assessment:

Reassessment appointment should be made every **3 years**.

Hearing Aid Option:

Battery Size:

Programs:

Volume Control: Please refer to specific user guide

TO SWITCH YOUR AID ON-CLOSE THE BATTERY COMPARTMENT

TO SWITCH YOUR AID OFF-OPEN BATTERY COMPARTMENT

Regular Maintenance:

Open fit- Tube and dome are one unit and should be replaced every 6 months. **DOMES MUST NOT BE REMOVED BY PATIENT.**

RITE- a new dome should be put on every time the old one is removed. Old domes should be disposed of immediately. Wax guards should be replaced if aid not functioning well or obvious wax is seen.

Tubing in ear moulds should be replaced every 6 months

Red marker indicates right aid. Blue marker indicates left aid.

Contact Us:

- **(028) 9056 1307**- To arrange an appointment

**If you need any replacement items:**

**Batteries**

**Tubes (Regular/ Open fit)**

**Domes**

**Wax guards**

**TEXT- 07568102636**

Please include your **Name and Date of Birth** to ensure we can fulfil the request

**Alternatively:**

**MAIL-** Audiology Department, Ulster Hospital, Dundonald. BT16 1RH  
-To arrange an appointment or request batteries

**EMAIL:** [audiology.uhd@setrust.hscni.net](mailto:audiology.uhd@setrust.hscni.net) –To arrange an appointment or request batteries

Other Useful Agencies:

**Sensory Support Services:** Lisburn (028) 9260 7746/ Downpatrick (028) 4461 6915/ Ards (028) 9151 0136/ Belfast (028) 9504 0200

**British Tinnitus Association:** Telephone 0800 018 0527 Website: [www.tinnitus.org.uk](http://www.tinnitus.org.uk)

**Hearing Link:** Telephone: 0300 111 1113 Website: [www.hearinglink.org/contactus](http://www.hearinglink.org/contactus)

**Action on Hearing Loss:** Telephone: 0808 808 0123 [www.actiononhearingloss.org.uk](http://www.actiononhearingloss.org.uk)

---

Your hearing aid remains the property of the NHS.