

Paediatric Hearing Aid Care Leaflet

Date of Fitting:

Hearing Aid Option:

Battery Size:

Volume Control: Please refer to specific user guide

To turn the hearing aid ON- Close the battery compartment

To turn the hearing aid OFF- Open the battery compartment

Regular Maintenance:

Open fit- The tube and dome is one unit and should be replaced every 6 months. **DOMES MUST NOT BE REMOVED BY PATIENT.**

RITE- **Domes are single use items. If you change the wax guard you must also replace a new dome.** Wax guards should be replaced if aid is not functioning well or ceases to work. Old domes should be disposed of immediately.

Tubing in ear moulds should be replaced every 6 months either by yourself or an appointment can be made at the department.

Your child will require regular ear moulds whilst they are still growing. If their hearing aid starts whistling or the mould seems visibly slack in the ear please contact the department and arrange an appointment for an impression to be taken.

Red marker indicates right aid. Blue marker indicates left aid.

Contact Us:

(028) 9056 1307- To arrange an appointment

If you need any replacement items:

Batteries

Tubes (Regular/ Open fit)

Domes

Wax guards

TEXT- 07568102636

Please include your **Name and Date of Birth** and put **PAEDIATRIC** in the subject field to ensure we can fulfil the request

Alternatively:

MAIL- Audiology Department, Ulster Hospital, Dundonald. BT16 1RH
-To arrange an appointment or request batteries

EMAIL: audiology.uhd@setrust.hscni.net –To arrange an appointment or request batteries or:

For queries/questions only you can email
paediatric.audiology@setrust.hscni.net.

Other Useful Agencies:

National Deaf Children Society (www.ncds.org.uk):
Royal National Institute for Deaf People (rnid.org.uk)

Your hearing aid remains the property of the NHS.